

Apaf-1 Antibody [5E11]


CATALOG NUMBER: PM-2016

Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse, Rat
Immunogen	A peptide corresponding to amino acids near the carboxy terminus of human Apaf1.
Conjugate	Unconjugated
Tested Applications	ELISA, IHC-P
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Apaf-1 Monoclonal Antibody is immunoaffinity chromatography purified IgG.
Clonality	Monoclonal
Isotype	IgG1,k
Physical State	Liquid
Buffer	Apaf-1 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Apaf-1 monoclonal antibody can be stored at -20°C, stable for one year.



Apaf-1 Antibody [5E11] 1

Immunohistochemistry of Apaf1 in K562 cells with Apaf1 antibody at 0.5 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).