



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Hemagglutinin Antibody [1E7D8]


CATALOG NUMBER: PM-4008

Specifications

Host Species	Mouse
Species Reactivity	Virus
Immunogen	A peptide corresponding to 13 amino acids in the middle of the Hemagglutinin protein.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Hemagglutinin Monoclonal Antibody is Protein A purified.
Clonality	Monoclonal
Physical State	Liquid
Buffer	Hemagglutinin Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Hemagglutinin monoclonal antibody can be stored at -20°C, stable for one year.



Hemagglutinin Antibody [1E7D8] 1

Western blot analysis of (A) 5 ng and (B) 25 ng of recombinant HA1 with Hemagglutinin antibody at 1 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).