

PUMA Antibody [10D4G7]


CATALOG NUMBER: PM-4529

Specifications

Host Species	Mouse
Species Reactivity	Human, Rat
Immunogen	Recombinant protein corresponding to amino acids 76 - 170 of human PUMA-alpha.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	PUMA Monoclonal Antibody is Protein A purified.
Clonality	Monoclonal
Isotype	IgG1
Physical State	Liquid
Buffer	PUMA Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	PUMA monoclonal antibody can be stored at -20°C, stable for one year.



PUMA Antibody [10D4G7] 1

Western blot analysis of PUMA expression in K562 cell lysate with PUMA antibody at (A) 2.5 and (B) 5 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).