



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## EndoG Antibody [7G1C10]


CATALOG NUMBER: PM-4577

### Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein corresponding to amino acids 51 - 140 of human EndoG.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	EndoG Monoclonal Antibody is Protein A purified.
Clonality	Monoclonal
Isotype	IgG1
Physical State	Liquid
Buffer	EndoG Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	EndoG monoclonal antibody can be stored at -20°C, stable for one year.



EndoG Antibody [7G1C10] 1

Western blot analysis of EndoG expression in HepG2 cell lysate with EndoG antibody at (A) 2.5 and (B) 5 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).