



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

MyD88 Antibody [2E9C2]


CATALOG NUMBER: PM-4673

Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse
Immunogen	Mouse monoclonal MyD88 antibody was raised against a recombinant protein corresponding to amino acids 176 to 280 of human MyD88.
Conjugate	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	MyD88 Monoclonal Antibody is Protein A purified.
Clonality	Monoclonal
Isotype	IgG2b
Physical State	Liquid
Buffer	MyD88 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	MyD88 monoclonal antibody can be stored at -20°C, stable for one year.



MyD88 Antibody [2E9C2] 1

Western blot analysis of MyD88 in EL4 whole cell lysate with MyD88 antibody at 2 µg/mL.



MyD88 Antibody [2E9C2] 2

Immunohistochemistry of MyD88 in human heart with MyD88 antibody at 2.5 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).