

## Hax1a Antibody [8F9G7]


CATALOG NUMBER: PM-4717

### Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse, Rat
Immunogen	Mouse monoclonal Hax1a antibody was raised against a 15 amino acid synthetic peptide near the amino terminus of human Hax1a.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	Hax1a Monoclonal Antibody is Protein G purified.
Clonality	Monoclonal
Isotype	IgG1,k
Physical State	Liquid
Buffer	Hax1a Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Hax1a monoclonal antibody can be stored at -20°C, stable for one year.



**Hax1a Antibody [8F9G7] 1**

Western blot analysis of Hax1a in human brain tissue lysate with Hax1a antibody at (A) 1 and (B) 2  $\mu$ g/mL.

**Disclaimer**

<b>Disclaimer</b>	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).