

Hax1a Antibody [9F3D11]


CATALOG NUMBER: PM-4718

Specifications

Host Species	Mouse
Species Reactivity	Human, Rat
Immunogen	Mouse monoclonal Hax1a antibody was raised against a 15 amino acid synthetic peptide near the amino terminus of human Hax1a.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Hax1a Monoclonal Antibody is Protein G purified.
Clonality	Monoclonal
Isotype	IgG3
Physical State	Liquid
Buffer	Hax1a Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Hax1a monoclonal antibody can be stored at -20°C, stable for one year.



Hax1a Antibody [9F3D11] 1

Western blot analysis of Hax1a in human brain tissue lysate with Hax1a antibody at (A) 1 and (B) 2 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).