

## IRAK Antibody [8F1A7]


CATALOG NUMBER: PM-5317

### Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse
Immunogen	Recombinant protein corresponding to amino acids 618 - 712 of human IRAK.
Conjugate	Unconjugated
Tested Applications	ELISA, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


### Properties

Purification	IRAK Monoclonal Antibody is Protein G purified.
Clonality	Monoclonal
Isotype	IgG1
Physical State	Liquid
Buffer	IRAK Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	IRAK monoclonal antibody can be stored at -20°C, stable for one year.



**IRAK Antibody [8F1A7] 1**

Western blot analysis of IRAK in mouse liver lysate with IRAK antibody at (A) 1 and (B) 2  $\mu\text{g/mL}$ .



**IRAK Antibody [8F1A7] 2**

Immunohistochemistry of IRAK in rat liver tissue with IRAK antibody at 5  $\mu\text{g/mL}$ .

**Disclaimer**

<b>Disclaimer</b>	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).