



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Seasonal H1N1 Neuraminidase Antibody [10C5E11]


CATALOG NUMBER: PM-5917

Specifications

Host Species	Mouse
Species Reactivity	Virus
Immunogen	Mouse monoclonal neuraminidase antibody was raised against a synthetic peptide from the seasonal influenza (H1N1) neuraminidase protein.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Seasonal H1N1 Neuraminidase Monoclonal Antibody is Protein A purified.
Clonality	Monoclonal
Isotype	IgG2b
Physical State	Liquid
Buffer	Seasonal H1N1 Neuraminidase Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Seasonal H1N1 Neuraminidase monoclonal antibody can be stored at -20°C, stable for one year.



Seasonal H1N1 Neuraminidase Antibody [10C5E11] 1

Seasonal influenza A N1 antibody (Cat. No. PM-5917) specifically recognizes seasonal (H1N1) N1, and does not cross-react with peptide corresponding to swine-origin influenza A (S-OIV, H1N1) N1 pept...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).