



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## GAPDH Antibody [12D3D5]


CATALOG NUMBER: PM-7569

### Specifications

Host Species	Mouse
Species Reactivity	Human, Mouse, Rat
Immunogen	GAPDH antibody was raised against a 16 amino acid peptide near the carboxy terminus of human GAPDH.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	GAPDH antibody is human, mouse, and rat reactive.
Predicted Molecular Weight	Predicted: 37 kDa Observed: 36 kDa

### Properties

Purification	GAPDH antibody is Protein A purified.
Clonality	Monoclonal
Isotype	IgG
Physical State	Liquid
Buffer	GAPDH antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	GAPDH antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



GAPDH Antibody [12D3D5] 1

Western blot analysis of GAPDH in multiple cell and tissue lysates with GAPDH antibody at 1 µg/ml. Lanes 1-23, 293, A431, A549, Daudi, HeLa, HepG2, Jurkat, K562, MOLT, 3T3, Raji, Ramos, U937, human...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).