

Vinculin Antibody [8B5]


CATALOG NUMBER: PM-7813

Specifications

Host Species	Mouse
Species Reactivity	Chicken, Human, Mouse, Rabbit, Rat
Immunogen	Vinculin antibody was raised against a 23 amino acid peptide near the amino terminus of human Vinculin.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	Vinculin antibody is human, mouse, rat, rabbit and chicken reactive. At least three isoforms of Vinculin are known to exist.
Predicted Molecular Weight	Predicted: 117 kDa Observed: 120 kDa

Properties

Purification	Vinculin antibody is purified via Protein A purification.
Clonality	Monoclonal
Isotype	IgG
Physical State	Liquid
Buffer	Vinculin Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Vinculin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



Vinculin Antibody [8B5] 1

Western blot analysis of Vinculin in 293, A431, A549, HeLa, HepG2, K562, 3T3, Raji, U937 cell lysate and human kidney, human breast, mouse brain, mouse lung, mouse spleen, rat heart, rat lung, rat ...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).