



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Avian Flu (H5N1) Protein Detection Set

CATALOG NUMBER: PSI-1808

Specifications

Species Reactivity	Human
Immunogen	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Antibodies are supplied as affinity chromatography purified IgG.
Clonality	Polyclonal
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	Antibody 1 mg/mL Peptide 200 µg/mL
Storage Conditions	H5N1 antibodies and peptides should be stored at -20°C, stable for one year.

Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).