



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## West Nile virus (WNV) Protein Detection Set

CATALOG NUMBER: PSI-1809

### Specifications

Species Reactivity	Human
Immunogen	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	Antibodies are supplied as affinity chromatography purified IgG.
Clonality	Polyclonal
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	Antibody 1 mg/mL Peptide 200 µg/mL
Storage Conditions	Stable at 4°C for three months, store at -20°C for up to one year.

### Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).