



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

# Complete Loading Control Detection Set Alpha-Tubulin, Beta-Actin, GAPDH

CATALOG NUMBER: PSI-1838

## Specifications

<b>Species Reactivity</b>	Human
<b>Immunogen</b>	Rabbit polyclonal antibodies were raised against peptides corresponding to amino acid sequences from each of the corresponding proteins.
<b>Tested Applications</b>	WB
<b>User Note</b>	Optimal dilutions for each application to be determined by the researcher.

## Properties

<b>Purification</b>	Antibodies are supplied as affinity chromatography purified IgG.
<b>Clonality</b>	Polyclonal
<b>Physical State</b>	Liquid
<b>Buffer</b>	PBS containing 0.02% sodium azide.
<b>Concentration</b>	Antibody 1 mg/mL
<b>Storage Conditions</b>	Stable at 4°C for three months, store at -20°C for up to one year.

## Disclaimer

<b>Disclaimer</b>	This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).