

TIGIT Antibody [7E5]


CATALOG NUMBER: RF16051

Specifications

Host Species	Mouse
Species Reactivity	Human
Immunogen	TIGIT antibody was raised against the extracellular domain of human TIGIT.
Conjugate	Unconjugated
Tested Applications	ELISA, Flow, ICC, IF, IHC-P
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	Predicted: 26 kDa


Properties

Purification	TIGIT Antibody is supplied as protein A purified IgG2b.
Clonality	Monoclonal
Isotype	IgG2b
Physical State	Liquid
Buffer	TIGIT Antibody is supplied in PBS containing 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/mL
Storage Conditions	TIGIT antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



TIGIT Antibody [7E5] 1


Immunocytochemistry of TIGIT in over expressing HEK293 cells using TIGIT antibody and control mouse IgG antibody (left corner box) at 1 µg/ml.



TIGIT Antibody [7E5] 2

Immunofluorescence of TIGIT in over expressing HEK293 cells using TIGIT Antibody at 1 µg/ml.


Green: TIGIT Antibody [7E5] (RF16051)
Blue: DAPI staining



TIGIT Antibody [7E5] 3

Immunofluorescence of TIGIT in human stomach carcinoma tissue using TIGIT Antibody at 5 µg/ml.

Green: TIGIT Antibody [7E5] (RF16051)
Blue: DAPI staining



TIGIT Antibody [7E5] 4

Immunohistochemistry of TIGIT in human stomach carcinoma tissue using TIGIT Antibody and control mouse IgG (corner box) at 2 µg/ml.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).