

## PDL1 ELISA Matched Antibody Pair (Risk Free)

CATALOG NUMBER: RFP-0306

### Specifications

Host Species	Llama, Mouse
Species Reactivity	Human
Immunogen	PD-L1 antibodies were raised against the extracellular domain of human PD-L1.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	Mouse Monoclonal: Antibodies supplied as protein A purified IgG1. Single Domain: Antibodies supplied as affinity chromatography purified via nickel column.
Clonality	Monoclonal
Isotype	IgG1 and sdab
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	Antibody 1 mg/mL
Storage Conditions	Mouse Monoclonal: Antibodies can be stored at 4°C for three months and -20°C, stable for up to one year.  Single Domain: Antibody should be stored in working aliquots at -20C, stable for up to one year.  As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Disclaimer

Disclaimer	For research use only. For additional information, visit ProSci's <a href="#">Terms and Conditions Page</a> . This product is for research use only.
------------	---