



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

TIGIT ELISA Matched Antibody Pair (Risk Free)


CATALOG NUMBER: RFP-0501

Specifications

Host Species	Mouse
Species Reactivity	Human
Immunogen	TIGIT antibody was raised against the extracellular domain of human TIGIT
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	Antibodies are supplied as protein A purified IgG1.
Clonality	Monoclonal
Isotype	IgG1
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	Antibody 1 mg/mL
Storage Conditions	Antibodies can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



TIGIT ELISA Matched Antibody Pair
(Risk Free) 1

A sandwich ELISA was performed using the anti-TIGIT mAb RF16055 (5 µg/ml) as the capture antibody. Biotin-labeled anti-TIGIT mAb RF16058-biotin (1 µg/ml) and streptavidin-HRP (0.1 µg/ml) were used...

Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).