

BAP3 Antibody


CATALOG NUMBER: 4503

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	BAP3 antibody was raised against a 12 amino acid synthetic peptide from near the amino terminus of human BAP3. The immunogen is located within the first 50 amino acids of BAP3.
Conjugate	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	BAP3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	BAP3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	BAP3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.




BAP3 Antibody 1

Western blot analysis of BAP3 in SK-N-SH cell lysate with BAP3 antibody at 2 µg/mL.



BAP3 Antibody 2

Immunohistochemical staining of human brain tissue using BAP3 antibody at 2.5 µg/mL.



BAP3 Antibody 3

Immunofluorescence of BAP3 in human brain tissue with BAP3 antibody at 20 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).