

TNPO3 Antibody


CATALOG NUMBER: 4599

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	TNPO3 antibody was raised against a 18 amino acid synthetic peptide from near the amino terminus of human TNPO3. The immunogen is located within amino acids 220 - 270 of TNPO3.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	Predicted: 95, 101 kDa Observed: 96 kDa


Properties

Purification	TNPO3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	TNPO3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	TNPO3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



TNPO3 Antibody 1

Western blot analysis of TNPO3 in rat liver tissue lysate with TNPO3 antibody at (A) 1 and (B) 2 µg/mL.



TNPO3 Antibody 2

Western blot analysis of TNPO3 expression in Jurkat cell lysate with TNPO3 antibody at (A) 1 and (B) 2 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).